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
### General instructions for use and safety recommendations for the application of dental instruments

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#### Area of application

These general instructions for use and safety recommendations apply to all products and have to be generally observed. Non-observance of these instructions for use and safety recommendations increases the risk of injury and may impair the proper function of the instruments.

Separate instructions for use are enclosed in the packaging of products that require more detailed information (packaging is provided with the  symbol). These take precedence over the general instructions.

#### Storage

Store wrapped devices protected from UV radiation and high temperatures in a clean and dry environment. Do not store in the same room as solvents or chemicals.

Medical devices without direct identification marking have to be stored with their original packaging.

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### 1. Proper use

- Those instruments that are supplied non-sterile have to be prepared prior to first use.
- Make sure that the power systems (hand-pieces and contra-angles) used are clean and perfect from a technical and hygienic aspect.
- Depending on their shank type, insert the instrument into the chuck as deeply as possible. Make sure that they are properly locked.
- The instrument must be rotating at the desired speed before contact is made with the tissue or material.
- Avoid jamming and using the instrument as a lever. Excessive contact pressure has to be avoided.
- For eye protection wear safety glasses. Use appropriate respiratory protection (mouth and nose). In the dental laboratory, use appropriate suction unit.
- Avoid contact with the instruments' working parts as this may increase the risk of injury.
- Make sure to provide sufficient cooling by means of air/water spray. Additional external cooling is required when using instruments with extra-long shanks or oversized working parts.
- Damaged or corroded instruments have to be discarded.

### 2. Recommended speeds

Please make sure to observe the recommendations for use and recommended speeds as indicated in the instructions for use and on the packaging of the products.

- $\varnothing_{max}$  300 000 rpm means: Suited for micro-motor contra-angles and turbines with stable ball bearings. Do not use in turbines with air bearing.
- $\varnothing_{max}$  200 000 rpm means: Suited for micro-motor hand pieces and contra-angles or lab hand pieces up to the speed indicated. Do not use in turbines.

Not observing the maximum permissible speed leads to an increased safety risk.

### 3. Contact pressure

Excessive contact pressure (> 2N) has to be avoided.

- In cutting instruments, this can lead to damage to the working part and to chipping of the blades as well as an excessive generation of heat.
- In abrasive instruments, increased contact pressure may lead to stripping of the grit or to clogging of the instruments and increased heat generation.
- Increased contact pressure may also lead to thermal damage to the pulp or, in case of damaged blades, to rough surfaces.
- In the extreme cases, instrument breakage may even occur.

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#### 4. Guideline on the number of times rotary instruments can be used

The below values are guidelines. The service life of instruments may differ from these values as this depends on the application and/or the material treated.

Stainless steel instruments:	up to 4 x
Tungsten carbide instruments:	up to 15 x
Diamond instruments:	up to 25 x
Polishers and ceramic abrasives:	up to 10 x
Endodontic instruments:	
Wide canals:	- max. 8 x
Average canals:	- max. 4 x
Narrow canals:	just use 1 x

The reuse of disposable products is not permitted. These are intended for single use and cannot be reprocessed.

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#### 5. Disposal

The devices have to be disposed of in tight, non-breakable and puncture-proof containers (protection from contamination).

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#### 6. Disinfection, cleaning and sterilisation

Those instruments that are supplied non-sterile have to be prepared prior to first use. Information on reprocessing, as amended, can be downloaded from our website <http://ifr.drendel.com> or ordered from the manufacturer Drendel + Zwilling DIAMANT GmbH.

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#### 7. Specific instructions for individual instrument types

- Avoid any contact with H<sub>2</sub>O<sub>2</sub> (hydrogen peroxide). Make sure that the specified immersion times in the cleaning and disinfecting agents are not exceeded. The carbide working parts would be attacked (a typical indication is black staining of the instrument) reducing the instrument's service life.
- Tool steel instruments corrode and can therefore not be sterilized with a sterilization method using saturated steam (e. g. autoclave).
- To achieve an optimal surface roughness, subsequent finishing is necessary after using a diamond grinding instrument with coarse or very coarse grit.
- Use a disc guard for diamond discs when working intraorally.
- When using diamond separating strips and diamond strips please avoid contact with the gingiva as there is a risk of injury.
- To avoid overstressing of the instrument, root canal instruments made of nickel titanium have to be used in a torque limited motor.
- Stainless steel root canal reamers (type "Gates", "Müller") are only intended for preparation of the coronal portion of the canal.
- Clean and sharpen Sintered Diamonds every now and then with the cleaning stone.
- Always use polishers and brushes with low pressure to minimize heat generation. Always polish in circular motion. To achieve brilliant high shine, use the polishers in the indicated sequence when using multiphase polishing systems.
- Polishers and Arkansas abrasives have to be prepared with suitable, alcohol free agents (e.g. DC EVO, Komet Dental/Alpro Medical).
- Surgical instruments with FG shank may only be used in a micromotor (red contra-angle). They are not authorized for use in a dental tubine.

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#### 8. Potential effects of using worn instruments

- Discard any instruments with damaged or worn working parts as the use of damaged or worn instruments would have a negative effect on the work result.
- Damaged and deformed cutting blades will cause the instrument to vibrate and lead to poor preparation margins and rough surfaces.
- Void spots on the surface of diamond instruments are a sign of missing diamond particles and reduce the instrument's abrasive efficiency. An inferior diamond coating quality will result in excessive heat generation. Excessive contact pressure or temperatures may cause damage to the tissue.

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## 9. Safety and possible side effects

The above mentioned recommendations with respect to cooling, contact pressure, disinfection, cleaning and sterilization are to be strictly observed. The instruments should only be used for the intended application, as per the symbolic identification. Non-observance of these safety recommendations may lead to damage of the power system and/or injury, such as thermal necrosis, undesired preparation of tissue, damage to tissue or nerves, violation of the biological width, or infections. During preparation, some instruments may generate metallic abrasion which may lead to the presence of artifacts in MRI diagnostics.

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## 10. Liability

The user is solely responsible for checking the devices prior to use to ensure that they are suited for the intended purposes. The correct application of the devices is the responsibility of the user and the operator. In case of contributory negligence by the user, Drendel + Zweiling DIAMANT GmbH partially or totally declines liability towards the user for all resulting damages, particularly if these are due to the non-observance of the recommendations for use or other misuse by the user. Store products out of children's reach. For dental use only.

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## 11. Contains hazardous substances

Devices marked with a CAS number (CAS: 7440-48-4) contain more than 0.1 % weight percent of cobalt. As a CMR substance of class IB, cobalt is classed as possibly carcinogenic, mutagenic and/or toxic to reproduction. Evidence has shown that the quantities of cobalt released by medical devices are so low that they do not pose a risk and no precautions have to be taken, provided that the medical device is used correctly and according to its intended purpose.

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## 12. Procedure in case of serious adverse events

In case of serious adverse events, please proceed in line with the legal regulations in your country.

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